



INSTRUCTIONS

Before use gently shake the bottle of reagent to clear any possible crystallization from the dripolator. A plastic pipette is included to assist if you have any dripolator blockage problems. To remove the dripolator carefully lever out with the edge of a knife. Make sure you wash the pipette thoroughly with water after each use to avoid any cross contamination of your reagents.

Mandelin reagent is a yellow/orange colour and may have a yellow precipitate, this is normal. Marquis and Mecke are clear and Simon's and Robatest have a light red tint.

For Marquis, Mandelin & Mecke .

1. Use a white ceramic plate.
2. Scrape a very small amount of powder from the edge of your pill onto the plate.
3. Remove the cap from the bottle of reagent and drip 1 or 2 drops onto the scraping.
4. Observe any colour change and then check the results chart to determine the possible substance in your pill.

For Simon's and Robatest.

1. Follow the above instructions and then add 2 drops of Buffer to the scraping as well, then check the results chart.
2. For Robatest use a larger scraping, around 10% of the pill for a more accurate result.

SAFETY

Make up a small bottle with about 10 grams of Sodium Bicarbonate (baking soda) dissolved in water and pour this on any reagent spills.

Marquis, Mandelin & Mecke reagents contain a concentrated acid, so please avoid contact with skin and eyes!!!